CE
0499

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CE
0499

DFL ALLOY
Dental Amalgam Alloy
English

SILVER FILINGS HIGH CONTENT OF COPPER

PROFESSIONAL USE ONLY.

DFL ALLOY contains: Silver (Ag) 45%, Tin (Sn) 31% and Copper (Cu) 24%.

DFL ALLOY is a silver fillings with high content of copper, superfine particles and excellent quality. It is made with pure metals. Their particles or micro granules are formed free from the contamination in the superficial layer, which results in a faster amalgamation with the mercury and consequently, larger mechanical resistance. The grind is made with easiness, without pressure and in an efficient way. The scientific mixture of the particle sizes in **DFL ALLOY**

facilitates the condensation done with small condensers and allows the production of denser restorations. The dust like micro particles (smaller than 8 micrometers) were eliminated, since they provoke a high flow and a great variation in the expansion during the binding. **DFL ALLOY** is virtually free from the Gamma II phase, responsible for the amalgam failures due to corrosion.

Indications:

DFL ALLOY is indicated for restorations in posterior teeth with little or any aesthetic compromising.

Presentations:

Bottle with 30 g of filings.
Pot with 50 capsules containing 320 mg of filings and 320 mg of mercury.
Pot with 50 capsules containing 400 mg of filings and 400 mg of mercury.
Pot with 50 capsules containing 500 mg of filings and 500 mg of mercury.
Pot with 50 capsules containing 600 mg of filings and 600 mg of mercury.
Pot with 500 capsules containing 320 mg of filings and 320 mg of mercury.
Pot with 500 capsules containing 400 mg of filings and 400 mg of mercury.
Pot with 500 capsules containing 500 mg of filings and 500 mg of mercury.
Pot with 500 capsules containing 600 mg of filings and 600 mg of mercury.

Storage and Conservation:

AMÁLGAMA DFL ALLOY must be stored in its original package and in a dry place. Keep the containers tightly closed after use. Under these conditions, the product is valid for up to 5 years after manufacturing date. Do not use this product if shelf life had expired.

Instructions for Use:

The use of gloves is recommended during the procedure.

Attention: The patient must be informed on the risk to occur sensitivity in the place of the procedure. This sensitization can be minimized through the use of proper isolation (rubber sheet) and pulp protection.

Dosage

As occurs for all the amalgams, the mercury-filings rate is very important. We recommended the rate of 1:1. For the use of filings powder, we recommend the use of mercury of good origin.

Manual grind

Put the mercury and the filings in a grail and triturate with a pistil for 50 seconds applying a pressure of approximately 1.20 Kg and 120 revolutions per minute. The amalgamation will be ready when a plastic, flat and shiny is obtained. To remove the excess of mercury, put the mixture on a clean and dry rubber sheet or on a cloth to squeeze. The removal of this excess should be made immediately before the insert of the amalgam in the cavity preparation

Mechanical grind :

Times suggested for **DFL ALLOY powder:**

Amalgamator adjusted for fixed or high-speed	Time
Silamat, Varimix	5 to 6 sec
Torit, Capmaster	8 to 10 sec
Dozamat, Dentomat, Aurimix	23 to 25 sec

Times suggested for **DFL ALLOY in capsules:**

Reference	Time	Speed
Varimix II	10 sec	High (H) 2 or 3
Varimix II	10 sec	High (H)
Ultramat 2 (SDI)	15 sec	-----
Capmix (Espe)	10 sec	-----
Amalga Mix (Gnatus)	20 sec	-----

The suggested times are approximate, once there is a variable in the production tolerance and in the electric currents, from place to place. As rule, the time of grind should be increased, when larger amounts of amalgam are prepared and reduced, when smaller amounts are desired.

Condensation and Sculpture

Absolute isolation should be used to avoid the contact of the amalgam with the oral cavity. Put the amalgam in the angles and in the retentions of the preparation, with pressure, using small condensers. Fill out the cavity with an excess, above the concave-superficial margin and remove the mercury that appears after the condensation and the burnishing, remove the matrix, sculpt the occlusal face, compose the margins and remove the cervical excess. If desired, an automatic condensation can be used.

Setting time: 8 to 12 minutes.

Finishing and Polishing

The final finishing and polishing should happen in a subsequent session. Use finishing multi-laminated drills to remove the marks of the instruments. Make the normal polishing using disks or rubber tips successively less abrasives, until obtaining the desired flat surface. It can be used pumice stone or polishing paste to obtain the final gloss, with the aid of felt tips. Product wastes must not be reused. They must be disposed as the local law in force.

Precautions:

The mercury leakage in the floor and in the pieces of furniture of the clinic should be avoided. The product should be always handled with gloves, mask and safety glasses. The clinic should have a high exhaustion system and failures in the high suction system should be avoided. Safe mechanical amalgamators, which do not present mercury leaks, should be used. Automatic condensers should not be used. The removal of restorations should be made with high suction, new drills and irrigation with cold water. In case of contact oral soft tissues, you should immediately remove the material with hydrophilic cotton soaked with water and rinse thoroughly with water. This product should not be used in pregnant women or in nursing mothers, due to the possibility of contamination of the fetus and of the breast milk by the mercury. This product must be used only by surgeon dentists and for the intended used described in the Instructions for Use. The empty capsules cannot be discarded in the environment, because they are polluted with mercury. They should be stocked in unbreakable and tightly closed containers, together with any amalgam residues generated by handled excess or by the sculpture of the material, and send, as soon as possible, to a chemical residues recovery laboratory or recycling.

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DFL Indústria e Comércio S.A.